


# How to collect stool (poop) using a bulk container

## 1 Label the container

- Your name
- Your date of birth
- Date & time of sample



## 2 Wash your hands



## 3 Urinate (pee) and flush




## 4 Collect stool using:

newspaper



flushable adhesive stool hat (included)



diaper



## 5 Put 2-3 spoonfuls of stool into container

Use a wooden stick (included) or a disposable spoon



Close lid tightly

## 6 Put container into biohazard bag – keep cool, not frozen



## 7 Throw out extra stool and supplies



## 8 Wash your hands



## 9 Call \_\_\_\_\_ at (\_\_\_\_)\_\_\_\_-\_\_\_\_\_ to return stool for testing